



Deep Clean and Standard Cleaning Procedure

After attending noninfectious patient					
Method	Standard				
New gloves should always be put on prior to cleaning. Remember to always wear gloves and other personal protection equipment, like a gown, cap or mask, as you could come in contact with blood or other body fluids during cleaning.	Follow correct glove use process				
Replace the used bed sheet immediately after each patient is dismissed. Put it in the clinical waste trashcans. Disinfect stretcher including rails and underside using disinfectant wipes.					
Disinfect the following using disinfectant wipes. Ensuring to wipe the surfaces evenly . Refer to CGP 215 Medical equipment cleaning/ disinfection for more details. -Medical equipment used such as thermometer, BP apparatus, pulse oximeter, LifePak, Stethoscope, etc. - Used MDT and Radios	Fully disinfected with concentrated level of approved disinfectant . The appropriate PPE must be used				
After disinfection, return items to designated place.					
Disinfect the exterior of BLS bag using disinfectant wipes. Return bag to ambulance cabinet.					
Disinfect using disinfectant wipes routinely touched areas such as door handles, etc.					
Dispose of any inventory items used or soiled as per policy.	Refer to OPP 120 Hazardous materials policy				
Dispose medical waste as per policy. Refill medical waste bag.	and procedures				
Remove gloves and mask. Hand hygiene.	Refer to CGP 129 Infection control program				
After attending a patient with suspected/ confirmed communicable dise	ase				
Method	Standard				
New gloves should always be put on prior to cleaning. Remember to always wear gloves and other personal protection equipment, like a gown, cap or mask, as you could come in contact with blood or other body fluids during cleaning.	Follow correct glove use process				
Prepare the cleaning items needed including mops. Prepare cleaning solution for mopping the ambulance floors.	For approved mop color and approved disinfectant, please refer to CGP 129 Infection control program				
Dispose of any inventory items used or soiled as per policy.	Refer to OPP 120 Hazardous materials policy and procedures				
Place stretcher in a clean, secured and shaded area (if possible) outside the ambulance. Ensure wheels are locked.	Keeping stretcher wheels locked ensure safety.				
Get the movable medical equipment from inside ambulance and place on top of stretcher	To facilitate cleaning inside ambulance				

Dispose of any inventory items used or soiled as per policy.	and procedures		
mopping the ambulance floors.	disinfectant, please refer to CGP 129 Infection control program Refer to OPP 120 Hazardous materials policy		
New gloves should always be put on prior to cleaning. Remember to always wear gloves and other personal protection equipment, like a gown, cap or mask, as you could come in contact with blood or other body fluids during cleaning. Prepare the cleaning items needed including mops. Prepare cleaning solution for	Follow correct glove use process For approved mop color and approved		
Method	Standard		
Patient's Cabin deep cleaning			
Once a month Deep Cleaning			
Update Operative IQ	For documentation		
attending a patient with suspected/ confirmed communicable disease			
Follow CGW 102 Dealing with spill workflow After dealing with spill, clean the ambulance following the procedures under "after	Refer to CGP 129 Infection control program		
Method	Standard		
When there is Body Fluid spill			
Update Operative IQ	For documentation		
Return stretcher inside the ambulance. Put new stretcher sheet. Hand hygiene.			
Remove gloves and mask. Hand hygiene.	Refer to CGP 129 Infection control program		
Dispose medical waste as per policy. Refill medical waste bag.	Refer to OPP 120 Hazardous materials policy and procedures		
Disinfect using disinfectant wipes routinely touched areas such as door handles, etc.			
Remove stretcher sheet. Disinfect stretcher including rails and underside using disinfectant wipes.			
Disinfect the exterior of BLS bag using disinfectant wipes. Return bag to ambulance cabinet.			
After disinfection, return items to designated place.	to wear appropriate PPEs.		
Refer to CGP 215 Medical equipment cleaning/ disinfection for more detailsMedical equipment used such as thermometer, BP apparatus, pulse oximeter, LifePak, Stethoscope, etc Used MDT and Radios	Other crew can do this simultaneously. Ensur		
Disinfect the following using disinfectant wipes. Ensuring to wipe the surfaces evenly.			
Mop ambulance patient treatment area floor using approved disinfectant with Blue mop. Allow to air dry for 10 minutes.	Follow recommended contact time of approved disinfectant solution used.		
 - Ambulance ceiling with AC grills and exhaust - Locked/ unopened lockers external area. For opened lockers, inside and outside. - Ambulance walls - Handles inside the ambulance - Seats - Fire extinguisher 	Fully disinfected with concentrated level of approved disinfectant . The appropriate PPE must be used.		
Disinfect the following using disinfectant wipes. Ensuring to wipe the surfaces evenly using appropriate number of wipes:			

	1		
Place stretcher in a clean, secured and shaded area (if possible) outside the ambulance. Ensure wheels are locked.	Keeping stretcher wheels locked ensure safety.		
Get the movable medical equipment from inside ambulance and place on top of stretcher	To facilitate cleaning inside ambulance		
Disinfect the following using National Ambulance approved disinfectant wipes. Ensuring to wipe the surfaces evenly using appropriate number of wipes: - Ambulance ceiling with AC grills and exhaust - Lockers inside and outside. - Ambulance walls - Handles inside the ambulance - Seats (seat belts and lap belt, under and behind - Fire extinguisher - Rear doors (wiring, small compartment and hinges)	Fully disinfected with concentrated level of approved disinfectant . The appropriate PPE must be used.		
Mop ambulance patient treatment area floor using approved disinfectant with Blue mop. Allow to air dry for 10 minutes.	Follow recommended contact time of approved disinfectant solution used.		
Disinfect the following using disinfectant wipes. Ensure to wipe all areas evenly. Refer to CGP 215 Medical equipment cleaning/ disinfection for more details. - Medical equipment such as thermometer, BP apparatus, pulse oximeter, LifePak, Stethoscope, Lucas, Entenox bag etc MDT - HX cylinder - Scoop, combi board and stair chair			
After disinfection, return items to designated place.	Other crew can do this simultaneously. Ensure		
Disinfect the exterior of BLS bag using disinfectant wipes. Return bag to ambulance cabinet.	to wear appropriate PPEs.		
Remove stretcher sheet. Disinfect stretcher including rails and underside using disinfectant wipes.			
Disinfect using disinfectant wipes routinely touched areas such as door handles, etc.			
Dispose waste as per policy. Disinfect trash bins. Refill waste bag.	Refer to OPP 120 Hazardous materials policy and procedures		
Remove gloves and mask. Hand hygiene.	Refer to CGP 129 Infection control program		
Return stretcher inside the ambulance. Put new stretcher sheet. Hand hygiene.	1		
Driver's cabin cleaning	let - de d		
Method	Standard		
Wear gloves, mask. Disinfect the following using National Ambulance approved disinfectant wipes. Ensuring to wipe the surfaces evenly:	_		
-Dashboard (cup holder, steering wheel, under the dashboard, compartment, abovethe sun visor, handle near door	Fully disinfected with concentrated level of approved disinfectant . The appropriate PPE		
 Middle compartment (cup holder), hand break Used MDT and Radios Wall behind seats and window Inside door and compartments on door 	must be used.		
Mop driver's cabin floor using the approved disinfectant with black/ gray mop. Allow to	Follow recommended contact time of approved		
air dry for 10 minutes.	disinfectant solution used.		
	Refer to CGP 129 Infection control program		
Remove gloves and mask. Hand hygiene.	incici to cor 123 infection control program		

Ambulance Exterior cleaning (minimum once per week)						
Method	Standard					
New gloves should always be put on prior to cleaning.	Follow correct glove use process					
Clean exterior including tyres, windshield, and step boards	Should be visibly clean with no blood or body					
elean exterior including tyres, windshield, and step sounds	substances					
First Aid Post Deep Clean (minimum monthly, or heavily soiled or Diseases known or suspected,)						
Method	Standard					
New gloves should always be put on prior to cleaning. Remember to always wear gloves and other personal protection equipment, like a gown, cap or mask, as you could come in contact with blood or other body fluids during cleaning.	Follow correct glove use process					
Replace the used bed sheeth immediately after each patient is dismissed. Put it in the clinical waste trashcans.						
Empty the trashcans, and wipe all of them with a disinfectant cleaner. Change the bags in every trashcan						
Wipe all the surfaces with an accepted disinfectant cleaner. This includes countertops, chairs, desk areas, patient tables, door handles and knobs, and any other surface that						
Clean all the equipment for patient care. This is medical tools, procedure equipment, blood pressure cuffs, thermometers, and any other portable equipment; these should be wiped with a disinfectant as they are used during the day.	Fully disinfected with concentrated level of approved disinfectant . The appropriate PPE must be used					
Wipe down of walls; interior cleaning of cabinetry and shelves.						
Tidy up the office waiting area. Make sure the trashcan is taken away and clean the chairs, table surfaces, and desk areas with an accepted disinfectant cleaner. Get rid of any garbage on the floor and vacuum or sweep the floors as much as needed.						
Remove gloves and apply sanitizer	Hand hygiene 99.9999% bacteria free					
Staff Rest Area, Medication Room (pharmacy) and Standby Stations whe	re we have full control - General Clean					
(minimum of once per week)						
Method	Standard					
New gloves should always be put on prior to cleaning	Follow correct glove use process					
Working top to bottom, disinfect key staff touch areas) eg sofa, table, using Alcohol wipes	Cleaned and staff touch points disinfected					
Sweep floor and then mop with Black/Gray Mop and Hot water (it is not expected that approved disinfectant is used here as this is for general cleanliness). The regular cleaner is expected to use chemical cleaner for floor. Mop handles shall be designated with the color available in the market excluding (Yellow, Red, Green and Blue)	Refer to CGP 129 Infection control program					
(Tellow, Neu, Green and Blue)						

Deep Clean and Standard Cleaning Procedure

Remove gloves and apply sanitizer

CGP153

Sep-22

Version 5.1





Mark on Cleaning Sign Sheet that the clean has been conducted

<u>п</u> п ш

Hand hygiene 99.9999% bacteria free